

MS3-MaxiSpin8

Mini Spin Centrifuge

Mini Spin Centrifuge is perfect designed with novel and unique appearance.

The instrument equipped with combined rotor and several types of tube work for 1.5ml/ 2.0ml / 0.5ml / 0.2ml tubes, 0.2ml x 8 PCR tube strips at the one cycle.

Mini spin centrifuge is designed user-friendly. Open/close the cover will automatically start/stop the centrifuge. Instrument built-in digital timing function.

Specification	MS3-MaxiSpin8
Speed range	6000 rpm
RCF	Approx. 2300 xg
Rotor capacity	8 x 1.5/2ml centrifuge tubes 8 x 0,5ml tubes 8 x 0,2ml tubes 2 x 8 x 0.2ml PCR strips
Digital timer set range	1-9999s or touch operation
Working noise	Maximum 45dB
Voltage	AC 220-230V, 50/60Hz
Power	45W
Fuse	250V, 1A , Ø5x20
Dimensions, WxDxH	156x176x121 mm
Net weight	1.5 kgs

