

MS2-MaxiFlake60

Laboratory Automatic Flake Ice Maker, Capacity: 60kg/day

Laboratory Automatic Flake Ice Maker is compact in size, made of highly durable material with energy efficient features. Increased efficiency gives quick production of ice flakes with ice making capacity 60 kg/day. Self contained storage bin eliminates the need of a separate storage unit.

Used in Pharmaceutical, Biochemical and biological studies, Biotechnology, Medical laboratories, Research laboratories and molecular biology laboratories.



- Produces granular shape ice flakes,
- Full automatic running and easy to use,
- Indicator system to indicate ice and water level,
- High efficiency CFC free compressor, economical and low energy consumption,
- Air cooled condenser,
- High quality stainless steel construction,
- Overload protection, low noise level,
- Ergonomically designed ice storage bin door,
- Product includes ice scoop, tap water hose and waste water hose.



Specifications

Model: MS2-MaxiFlake60

Flake ice production capacity	60kg/day
Flake ice storage capacity	15kg
Dimension (W x D x H)	398 x 542 x 700mm
Weight	40kg
Refrigerant gas	R290
Electrical voltage	220V/ 50Hz