


## MS3-MaxiPlate2 Microplate Centrifuge

Microplate Centrifuge is specially designed for the 96-hole microplates. It can quickly spin down droplets and condensation in the microplate.

Microplate centrifuge instrument have suitable for skirted, semi-skirted, non-skirted and all standard PCR plates.

The instrument have strong centrifugal power. Accelerate to 2500 rpm in a short time and from top speed to stop time around 10 seconds.

Smart size, compact operation. The cover can be fully opened to 95°, convenient to put and fetch the microplates.



Specification	MS3-MaxiPlate2
Speed range	2500 rpm
RCF	627 xg
Sample capacity	2 x PCR plates
Rotor	Fixed vertically
Digital timer set range	1-99s or touch operation
Voltage	DC 24V power adapter (AC 100V-230V, 50/60Hz)
Power	40W
Dimensions, WxDxH	190x219x186 mm
Net weight	2.7 kgs